

Product Datasheet

Recombinant human MUC16 Protein, C-His-Avi (HEK293) (orb1516600)

Catalog Number	orb1516600
Category	Proteins
Description	Recombinant human MUC16 Protein, C-His-Avi (HEK293), with a concentration of >0.5mg/ml, expressed in HEK293 cell, with a molecular weight of 40-70 kDa, and a purity of > 90% as determined by SDS-PAGE. It was purified by AC, and has a C-His-Avi tag.
Target	Mucin-16
Tag	C-His-Avi
Form/Appearance	Lyophilized or Liquid
Concentration	>0.5mg/ml
Buffer/Preservatives	PBS (pH7.4).
Reconstitution	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Purity	> 90% as determined by SDS-PAGE
Purification	AC
Protein Sequence	12660-12923/14507
MW	40-70 kDa
Expression System	HEK293 cell
Endotoxins	Not analyzed

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Note For research use only

Expiration Date 6 months from date of receipt.



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)