


Product Datasheet

Recombinant human IL-6 protein, C-His (HEK293) (orb1516624)

Catalog Number	orb1516624
Description	Recombinant human IL-6 protein, C-His (HEK293)
Target	Interleukin-6
Form/Appearance	Liquid
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	30-212/212
Purity	> 90% as determined by SDS-PAGE method
MW	20.8 kDa
Expiration Date	6 months from date of receipt.

Human IL-6R Bind with Human IL-6



Measured by its binding ability in a direct ELISA. Immobilized Human IL-6 (orb1516624) at 2 µg/mL (100 µl/well) can bind Human IL-6R-HRP, the minimum detection concentration is 7.81 ng/ml.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)