

Product Datasheet

Human Prostate Specific Antigen protein (orb1516687)

Catalog Number	orb1516687
Category	Proteins
Description	Human Prostate Specific Antigen protein, with a molecular weight of 34 kDa, and a purity of >95% as determined by SDS-PAGE.
Target	Prostate-specific antigen
Form/Appearance	Lyophilized or Liquid
Purity	>95% as determined by SDS-PAGE
Protein Sequence	Purified native protein
MW	34 kDa
Endotoxins	Not analyzed
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)