

Product Datasheet

Recombinant human Bax protein, DsbC & His (orb1516786)

Catalog Number	orb1516786
Category	Proteins
Description	Recombinant human Bax protein, DsbC & His, with a concentration of >0.5 mg/ml, expressed in E.coli, with a molecular weight of 49.5 kDa, and a purity of >95% as determined by SDS-PAGE. It was purified by AC, and has a DsbC & His tag.
Target	Apoptosis regulator BAX
Tag	DsbC & His
Form/Appearance	Lyophilized or Liquid
Concentration	>0.5 mg/ml
Buffer/Preservatives	20mM Tris-HCl (pH8.0).
Reconstitution	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Purity	>95% as determined by SDS-PAGE
Purification	AC
Protein Sequence	1-192/192
MW	49.5 kDa
Expression System	E.coli
Endotoxins	Not analyzed

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage


Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Note

For research use only

Expiration Date

6 months from date of receipt.



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)