


Product Datasheet

Bovine Fibrinogen (from plasma) (orb1516930)

Catalog Number	orb1516930
Category	Proteins
Description	Bovine Fibrinogen (from plasma), with a molecular weight of 340 kDa.
Target	Fibrinogen alpha chain
Form/Appearance	Lyophilized
Protein Sequence	Purified native protein
MW	340 kDa
Endotoxins	Not analyzed
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558