

Product Datasheet

Goat anti-Human C-reactive protein antibody (HRP) Antibody (orb1517273)

Catalog Number	orb1517273
Description	Goat polyclonal antibody against C-reactive protein conjugated to HRP
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Human C-reactive protein antibody (HRP)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Concentration	0.2mg in 0.5mg/ml
Note	For research use only
Application notes	WB: Approx. 25kDa band observed in Human Liver and Serum Lysates. Recommended concentration: 0.3µg/ml. Primary incubation was 1 hour.
Protein Sequence	Purified Human C reactive protein
MW	25 kDa
Entrez	1401
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)