

Product Datasheet

Pig Albumin ELISA Kit (orb1518893)

Catalog Number	orb1518893
Category	Assays and Kits
Description	Pig Albumin ELISA Kit
Target	Pig Albumin
Reactivity	Porcine
UniProt ID	P08835
Tested applications	ELISA
Dilution range	WB: , IP: , IHC: , ICC: , IHC-IF: , ICC-IF: , ChIP: , ChIP-chip: , ChIP-PCR: , ChIP-Seq: , ELISA: , IEP: , DD: , IF:
Assay Type	Sandwich
Range	0.69 - 500 ng/ml
Sample Types	serum_plasma_urine
Storage	2 - 8°C
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)