

Product Datasheet

Rabbit anti-CD31 Recombinant Monoclonal Antibody (orb1519538)

Catalog Number	orb1519538
Category	Antibodies
Description	CD31/PECAM-1 Antibody
Target	CD31/PECAM-1
Clonality	Recombinant
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Mouse
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Buffer/Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.09% Sodium Azide, rAlbumin-Free
Immunogen	Between 1 and 50
UniProt ID	Q08481
Tested applications	FC, IHC, WB
Dilution range	WB - 1:1000; IHC - 1:100 to 1:500. Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE tissue sections.
Application notes	Format: Whole IgG

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom


Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Antibody Type	Primary Antibody
Clone Number	BLR180J
Storage	2 - 8°C
Note	For research use only
NCBI	NP_001292086.1



Detection of mouse CD31 (shaded) in Ba/F3 cells by flow cytometry. Antibody: Rabbit anti-CD31 recombinant monoclonal antibody (orb1519538) or isotype control (unshaded). Secondary: DyLight® 650-conjugated goat anti-rabbit IgG.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)