

Product Datasheet

Donkey anti-Goat IgG Heavy and Light Chain Cross-Adsorbed Antibody (Biotin) (orb1520039)

Catalog Number	orb1520039
Category	Antibodies
Description	Goat IgG-Heavy and Light Chain Antibody
Target	Goat IgG-Heavy and Light Chain
Clonality	Polyclonal
Species/Host	Donkey
Conjugation	Biotin
Reactivity	Goat
Predicted Reactivity	Gallus, Human, Mouse, Porcine, Rabbit, Rat
Form/Appearance	Whole IgG
Concentration	0.5 mg/ml
Buffer/Preservatives	Phosphate Buffered Saline (PBS) containing 0.09% Sodium Azide
Tested applications	ELISA, ICC, IHC, WB
Dilution range	WB - 1:10,000 - 1:200,000; IHC - 1:250 - 1:2,500; ICC - 1:100 - 1:500; ELISA - 1:10,000 - 1:200,000
Application notes	Format: Whole IgG
Antibody Type	Secondary Antibody
Storage	2 - 8°C

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)