

Product Datasheet

Rabbit anti-Goat IgG Heavy and Light Chain Cross-Adsorbed Antibody (HRP) (orb1520042)

Catalog Number	orb1520042
Category	Antibodies
Description	Goat IgG-Heavy and Light Chain Antibody
Target	Goat IgG-Heavy and Light Chain
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	HRP
Reactivity	Goat
Predicted Reactivity	Equine, Gallus, Human, Mouse, Porcine, Rat
Form/Appearance	Whole IgG
Concentration	0.5 mg/ml
Buffer/Preservatives	Phosphate Buffered Saline (PBS) containing 0.05% Pro-Clean 400
Tested applications	ELISA, ICC, IHC, WB
Dilution range	WB - 1:5,000 - 1:50,000; IHC - 1:200 - 1:5,000; ICC - 1:200 - 1:5,000; ELISA - 1:10,000 - 1:100,000
Application notes	Format: Whole IgG
Antibody Type	Secondary Antibody
Storage	2 - 8°C

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)