

Product Datasheet

Piccolo Antibody (FITC) (orb152169)

Catalog Number	orb152169
Category	Antibodies
Description	Rabbit polyclonal to Picollo (FITC). Piccolo, also referred to as Aczonin, is a large protein which consists of an N-terminal Zn²⁺ finger, several piccolo-bassoon homology domains and C-terminal PDZ and C2 domains. In general it is found together with bassoon, a related huge multi-domain protein of the CAZ (cytoskeletal matrix assembled at active zones). Piccolo is a scaffolding protein for proteins involved in endo- and exocytosis of synaptic vesicles. Piccolo has also been shown to interfere with clathrin mediated endocytosis by binding to the F-actin and dynamin binding protein Abp1..
Target	Piccolo
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	FITC
Reactivity	Human, Mouse, Rat
Concentration	1 mg/ml
Buffer/Preservatives	640.91mM DMSO, 136.36 mM Ethanolamine, 126.89 mM chlorides, 9.09mM phosphates, 9.09mM NaHCO ₃
Purification	Protein A Purified
Immunogen	Recombinant Full length Protein, BC001304.1
UniProt ID	Q9Y6V0
MW	550kDa

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom


Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC


68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Tested applications	IHC, WB
Dilution range	WB (1:1000), IHC (1:400)
Application notes	1 µg/ml of SPC-197 was sufficient for detection of Piccolo in 10 µg of rat brain tissue lysate by colorimetric immunoblot analysis using goat ant rabbit IgG:HRP as the secondary antibody.
Specificity	Detects ~550kDa. No cross-reactivity to Bassoon.
Storage	Conjugated antibodies should be stored according to the product label
Note	For research use only
Entrez	27445
NCBI	NP_055325.2



Immunohistochemistry analysis using Rabbit Anti-Piccolo Polyclonal Antibody. Tissue: Muscle. Species: Mouse. Primary Antibody: Rabbit Anti-Piccolo Polyclonal Antibody at 1:400. Secondary Antibody: Alexa Fluor 488 Goat Anti-Rabbit. Counterstain: BTX (red). Localization: Selective staining of the NMJ.



Western blot analysis of Mouse, Rat brain cell lysates showing detection of Piccolo protein using Rabbit Anti-Piccolo Polyclonal Antibody. Primary Antibody: Rabbit Anti-Piccolo Polyclonal Antibody at 1:1000.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558