

Product Datasheet

Rabbit anti-DCAMKL1 Antibody (orb1523034)

Catalog Number	orb1523034
Category	Antibodies
Description	DCAMKL1 Antibody
Target	DCAMKL1
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human, Mouse
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Buffer/Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Immunogen	Between 150 and 200
UniProt ID	O15075
Tested applications	IP, WB
Dilution range	WB - 1:500 - 1:2,500; IP - 2 - 10 µg/mg lysate
Application notes	Format: Whole IgG
Antibody Type	Primary Antibody
Storage	2 - 8°C

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Note For research use only

NCBI **O15075.2**



Detection of human and mouse DCAMKL1 by western blot. Samples: Whole cell lysate (50 µg) prepared using NETN buffer from HeLa, HEK293T, Jurkat, and mouse TCMK-1 cells. Antibodies: Affinity purified rabbit anti-DCAMKL1 antibody orb1523034 used for WB at 1 µg/ml. Detection: Chemiluminescence with an exposure time of 3 minutes.



Detection of human DCAMKL1 by western blot of immunoprecipitates. Samples: Whole cell lysate (1 mg for IP; 20% of IP loaded) prepared using NETN buffer from 293T and cells. Antibodies: Affinity purified rabbit anti-DCAMKL1 antibody orb1523034 used for IP at 6 µg per reaction.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558