

Product Datasheet

Rabbit anti-RBM22 Antibody (orb1523925)

Catalog Number	orb1523925
Category	Antibodies
Description	RBM22 Antibody
Target	RBM22
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human, Mouse
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Buffer/Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Immunogen	Between 370 and C-terminus
UniProt ID	Q9NW64
Tested applications	IP, WB
Dilution range	WB - 1:2,000 - 1:10,000; IP - 2 - 10 µg/mg lysate
Application notes	Format: Whole IgG
Antibody Type	Primary Antibody
Storage	2 - 8°C

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Note For research use only

NCBI [NP_060517.1](#)



Detection of human and mouse RBM22 by western blot. Samples: Whole cell lysate from HEK293T (15 and 50 μ g) , HeLa (50 μ g) , Jurkat (50 μ g) , and mouse NIH 3T3 (50 μ g) cells. Antibodies: Affinity purified rabbit anti-RBM22 antibody orb1523925 used for WB at 0.1 μ g/ml. Detection: Chemiluminescence with an exposure time of 30 seconds.



Detection of human RBM22 by western blot of immunoprecipitates. Samples: Whole cell lysate (1 mg for IP; 20% of IP loaded) from HEK293T cells. Antibodies: Affinity purified rabbit anti-RBM22 antibody orb1523925 used for IP at 6 μ g/mg lysate.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558