

Product Datasheet

Rabbit anti-RBM22 Antibody (orb1523927)

Catalog Number	orb1523927
Category	Antibodies
Description	RBM22 Antibody
Target	RBM22
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human, Mouse
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Buffer/Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Immunogen	Between 175 and 225
UniProt ID	Q9NW64
Tested applications	IP, WB
Dilution range	WB - 1:10,000 - 1:25,000; IP - 2 - 10 µg/mg lysate
Application notes	Format: Whole IgG
Antibody Type	Primary Antibody
Storage	2 - 8°C

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Note For research use only


NCBI [NP_060517.1](#)



Detection of mouse RBM22 by western blot. Samples: Whole cell lysate (10 µg) from NIH 3T3, CT26, CH27, TCMK-1, and BW5147.3 cells prepared using NETN lysis buffer. Antibody: Affinity purified rabbit anti-RBM22 antibody (orb1523927) used for WB at 0.04 µg/ml. Detection: Chemiluminescence with an exposure time of 1 second.



Detection of human RBM22 by western blot. Samples: Whole cell lysate (10 µg) from HEK293T, HeLa, Jurkat, OVCAR-8, and SR cells prepared using NETN lysis buffer. Antibody: Affinity purified rabbit anti-RBM22 antibody (orb1523927) used for WB at 0.04 µg/ml. Detection: Chemiluminescence with an exposure time of 3 seconds.



Detection of human RBM22 by western blot of immunoprecipitates. Samples: Whole cell lysate (1.0 mg per IP reaction; 5% of IP loaded) from HEK293T cells prepared using NETN lysis buffer. Antibodies: Affinity purified rabbit anti-RBM22 antibody (orb1523927) used for IP at 6 µg per reaction. RBM22 was also immunoprecipitated by a previous lot of this antibody (orb1523927) and a second antibody against a different epitope of RBM22.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558