

## Product Datasheet

### Rabbit anti-Filamin B Antibody (orb1527905)

|                     |                                                                      |
|---------------------|----------------------------------------------------------------------|
| Catalog Number      | orb1527905                                                           |
| Description         | Filamin B Antibody                                                   |
| Species/Host        | Rabbit                                                               |
| Reactivity          | Human                                                                |
| Conjugation         | Unconjugated                                                         |
| Tested Applications | WB                                                                   |
| Immunogen           | Between 1250 and 1300                                                |
| Target              | Filamin B                                                            |
| Preservatives       | Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide |
| Form/Appearance     | Whole IgG                                                            |
| Concentration       | 200 µg/ml                                                            |
| Storage             | 2 - 8°C                                                              |
| Note                | For research use only                                                |
| Application notes   | Format: Whole IgG                                                    |
| Clonality           | Polyclonal                                                           |
| Antibody Type       | Primary Antibody                                                     |
| Uniprot ID          | <b>O75369</b>                                                        |
| NCBI                | <b>NP_001448.2</b>                                                   |
| Dilution Range      | WB - 1:2,000 - 1:10,000; IP - Not recommended                        |

#### Biorbyt Ltd.


7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

**Expiration Date**

12 months from date of receipt.



Detection of human Filamin B by western blot. Samples: Whole cell lysate from HeLa (5, 15 and 50  $\mu$ g) and HEK293T (T; 50  $\mu$ g) cells. Antibodies: Affinity purified rabbit anti-Filamin B antibody orb1527905 used for WB at 0.04  $\mu$ g/ml. Detection: Chemiluminescence with an exposure time of 30 seconds.

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)