

## Product Datasheet

### Rabbit anti-BRIP1/BACH1 Antibody (orb1529990)

|                     |                                                                        |
|---------------------|------------------------------------------------------------------------|
| Catalog Number      | orb1529990                                                             |
| Description         | BRIP1 Antibody                                                         |
| Species/Host        | Rabbit                                                                 |
| Reactivity          | Human                                                                  |
| Conjugation         | Unconjugated                                                           |
| Tested Applications | IP                                                                     |
| Immunogen           | Between 1125 and 1175                                                  |
| Target              | BRIP1                                                                  |
| Preservatives       | Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide |
| Form/Appearance     | Whole IgG                                                              |
| Concentration       | 1000 µg/ml                                                             |
| Storage             | 2 - 8°C                                                                |
| Note                | For research use only                                                  |
| Application notes   | Format: Whole IgG                                                      |
| Clonality           | Polyclonal                                                             |
| Antibody Type       | Primary Antibody                                                       |
| Uniprot ID          | <b>Q9BX63</b>                                                          |
| NCBI                | <b>NP_114432.1</b>                                                     |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

#### Biorbyt LLC.


68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

**Dilution Range**

WB - Not recommended. Use rabbit anti-BRIP1(BACH1) antibody A300-561A.; IP - 1 - 4 µg/mg lysate

**Expiration Date**

12 months from date of receipt.



Detection of human BRIP1/BACH1 by western blot of immunoprecipitates. Samples: Whole cell lysate from HeLa cells. Antibody: Affinity purified rabbit anti-BRIP1/BACH1 antibody.

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)