

Product Datasheet

Goat anti-Donkey IgG Heavy and Light Chain Antibody (Biotin) (orb1530946)

Catalog Number	orb1530946
Category	Antibodies
Description	Donkey IgG-Heavy and Light Chain Antibody
Target	Donkey IgG-Heavy and Light Chain
Clonality	Polyclonal
Species/Host	Goat
Conjugation	Biotin
Reactivity	Donkey
Form/Appearance	Whole IgG
Concentration	1 mg/ml
Buffer/Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.09% Sodium Azide
Tested applications	ELISA, ICC, IHC, WB
Dilution range	WB - 1:20,000 - 1:400,000; IHC - 1:500 - 1:5,000; ICC - 1:200 - 1:1,000; ELISA - 1:20,000 - 1:400,000
Application notes	Format: Whole IgG
Antibody Type	Secondary Antibody
Storage	2 - 8°C
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Expiration Date

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)