

Product Datasheet

Rabbit anti-Donkey IgG Heavy and Light Chain Antibody (Biotin) (orb1530970)

Catalog Number	orb1530970
Category	Antibodies
Description	Donkey IgG-Heavy and Light Chain Antibody
Target	Donkey IgG-Heavy and Light Chain
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Biotin
Reactivity	Donkey
Form/Appearance	Whole IgG
Concentration	1 mg/ml
Buffer/Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.09% Sodium Azide
Tested applications	ELISA, ICC, IHC, WB
Dilution range	WB - 1:20,000 - 1:400,000; IHC - 1:500 - 1:5,000; ICC - 1:200 - 1:1,000; ELISA - 1:20,000 - 1:400,000
Application notes	Format: Whole IgG
Antibody Type	Secondary Antibody
Storage	2 - 8°C

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Note

For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)