

Product Datasheet

Goat anti-Rabbit IgG Heavy and Light Chain Highly Cross-Adsorbed Antibody (FITC) (orb1531013)

Catalog Number	orb1531013
Category	Antibodies
Description	Rabbit IgG-Heavy and Light Chain Antibody
Target	Rabbit IgG-Heavy and Light Chain
Clonality	Polyclonal
Species/Host	Goat
Conjugation	FITC
Reactivity	Rabbit
Predicted Reactivity	Human, Mouse
Form/Appearance	Whole IgG
Concentration	1 mg/ml
Buffer/Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.09% Sodium Azide
Tested applications	FC, ICC, IF, IHC
Dilution range	IHC - 1:100 - 1:500. For use with Formalin-Fixed Paraffin-Embedded Tissues; ICC - 1:100 - 1:500; IF - 1:100 - 1:500. For use with cells
Application notes	Format: Whole IgG
Antibody Type	Secondary Antibody
Storage	2 - 8°C

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)