

Product Datasheet

F(ab')₂ Donkey Rabbit IgG Heavy and Light Chain Cross-Adsorbed Antibody (orb1531025)

Catalog Number	orb1531025
Category	Antibodies
Description	Rabbit IgG-Heavy and Light Chain Antibody
Target	Rabbit IgG Heavy and Light Chain
Clonality	Polyclonal
Species/Host	Donkey
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Rabbit
Predicted Reactivity	Bovine, Gallus, Goat, Human, Mouse, Rat, Sheep
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Buffer/Preservatives	Phosphate Buffered Saline (PBS) containing 0.09% Sodium Azide
Purification	Antigen Affinity Purified
Tested applications	ELISA, ICC, IHC, WB
Dilution range	WB - 1:1,000 - 1:20,000; IHC - 1:200 - 1:2,000; ICC - 1:200 - 1:2,000; ELISA - 1:1,000 - 1:20,000; for coating plates 1:50 - 1:250

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Application notes	Format: F(ab') ₂
Antibody Type	Secondary Antibody
Storage	2 - 8°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)