

Product Datasheet

F(ab')₂ Goat anti-Rabbit IgG-F(ab')₂ Fragment Cross-Adsorbed Antibody (orb1531030)

Description	Rabbit IgG-F(ab) ₂ Fragment Antibody
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, WB
Target	Rabbit IgG-F(ab) ₂ Fragment
Preservatives	Phosphate Buffered Saline (PBS) containing 0.09% Sodium Azide
Form/Appearance	F(ab') ₂
Concentration	0.5 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: F(ab') ₂
Clonality	Polyclonal
Antibody Type	Secondary Antibody
Dilution Range	WB - 1:1,000 - 1:20,000; IHC - 1:200 - 1:2,000; ICC - 1:200 - 1:2,000; ELISA - 1:1,000 - 1:20,000; for coating plates 1:50 - 1:250
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)