

Product Datasheet

Goat anti-Rabbit IgG Heavy and Light Chain Cross-Adsorbed Antibody (HRP) (orb1531045)

Catalog Number	orb1531045
Category	Antibodies
Description	Rabbit IgG-Heavy and Light Chain Antibody
Target	Rabbit IgG-Heavy and Light Chain
Clonality	Polyclonal
Species/Host	Goat
Conjugation	HRP
Reactivity	Rabbit
Predicted Reactivity	Bovine, Equine, Gallus, Human, Mouse, Porcine, Rat
Form/Appearance	Whole IgG
Concentration	0.5 mg/ml
Buffer/Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.05% Pro-Clean 400
Tested applications	ELISA, ICC, IHC, WB
Dilution range	WB - 1:5,000 - 1:50,000; IHC - 1:200 - 1:5,000; ICC - 1:200 - 1:5,000; ELISA - 1:10,000 - 1:100,000
Application notes	Format: Whole IgG
Antibody Type	Secondary Antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage 2 - 8°C

Note For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)