

Product Datasheet

Goat Rabbit IgG Heavy and Light Chain Cross-Adsorbed Antibody (FITC) (orb1531046)

Catalog Number	orb1531046
Category	Antibodies
Description	Rabbit IgG-Heavy and Light Chain Antibody
Target	Rabbit IgG Heavy and Light Chain
Clonality	Polyclonal
Species/Host	Goat
Isotype	IgG
Conjugation	FITC
Reactivity	Rabbit
Predicted Reactivity	Bovine, Equine, Gallus, Human, Mouse, Porcine, Rat
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Buffer/Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.09% Sodium Azide
Purification	Antigen Affinity Purified
Tested applications	FC, ICC, IF, IHC
Dilution range	IHC - 1:50 - 1:500; ICC - 1:50 - 1:500; IF - 1:50 - 1:500

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Application notes	Format: Whole IgG
Antibody Type	Secondary Antibody
Storage	2 - 8°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)