

Product Datasheet

Goat anti-Rabbit IgG Heavy and Light Chain Cross-Adsorbed Antibody (Cy5) (orb1531048)

Description	Rabbit IgG-Heavy and Light Chain Antibody
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Cy5
Tested Applications	FC, ICC, IF, IHC, WB
Target	Rabbit IgG-Heavy and Light Chain
Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	0.5 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Antibody Type	Secondary Antibody
Dilution Range	WB - 1:1,000 - 1:15,000. 5% non-fat dry milk in PBST or TBST is recommended for blocking and incubation of antibodies.; IHC - 1:50 - 1:500; ICC - 1:50 - 1:500; IF - 1:50 - 1:500
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)