

Product Datasheet

F(ab')₂ Goat anti-Rabbit IgG-F(ab')₂ Fragment Antibody (orb1531077)

Catalog Number	orb1531077
Category	Antibodies
Description	Rabbit IgG-F(ab) ₂ Fragment Antibody
Target	Rabbit IgG-F(ab) ₂ Fragment
Clonality	Polyclonal
Species/Host	Goat
Conjugation	Unconjugated
Reactivity	Rabbit
Form/Appearance	F(ab') ₂
Concentration	0.5 mg/ml
Buffer/Preservatives	Phosphate Buffered Saline (PBS) containing 0.09% Sodium Azide
Tested applications	ELISA, ICC, IHC, WB
Dilution range	WB - 1:1,000 - 1:30,000; IHC - 1:200 - 1:2,000; ICC - 1:200 - 1:2,000; ELISA - 1:1,000 - 1:30,000; for coating plates 1:100 - 1:500
Application notes	Format: F(ab') ₂
Antibody Type	Secondary Antibody
Storage	2 - 8°C

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Note

For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)