

Product Datasheet

Goat anti-Rabbit IgG-Fc Fragment Antibody (Biotin) (orb1531080)

Description	Rabbit IgG-Fc Fragment Antibody
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC, WB
Target	Rabbit IgG-Fc Fragment
Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Antibody Type	Secondary Antibody
Dilution Range	WB - 1:20,000 - 1:400,000; IHC - 1:500 - 1:5,000; ICC - 1:200 - 1:1,000; ELISA - 1:20,000 - 1:400,000
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)