

## Product Datasheet

### Goat anti-Rabbit IgM Heavy Chain Antibody (Biotin) (orb1531086)

|                            |                                                                                                     |
|----------------------------|-----------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit IgM Antibody                                                                                 |
| <b>Species/Host</b>        | Goat                                                                                                |
| <b>Reactivity</b>          | Rabbit                                                                                              |
| <b>Conjugation</b>         | Biotin                                                                                              |
| <b>Tested Applications</b> | ELISA, ICC, IHC, WB                                                                                 |
| <b>Target</b>              | Rabbit IgM-Heavy Chain                                                                              |
| <b>Preservatives</b>       | Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.09% Sodium Azide                     |
| <b>Form/Appearance</b>     | Whole IgG                                                                                           |
| <b>Concentration</b>       | 1 mg/ml                                                                                             |
| <b>Storage</b>             | 2 - 8°C                                                                                             |
| <b>Note</b>                | For research use only                                                                               |
| <b>Application notes</b>   | Format: Whole IgG                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                                          |
| <b>Antibody Type</b>       | Secondary Antibody                                                                                  |
| <b>Dilution Range</b>      | WB - 1:10,000 - 1:200,000; IHC - 1:250 - 1:2,500; ICC - 1:100 - 1:500; ELISA - 1:10,000 - 1:200,000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                     |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)