

Product Datasheet

Goat anti-Rat IgG Heavy and Light Chain Cross-Adsorbed Antibody (Cy5) (orb1531132)

Catalog Number	orb1531132
Category	Antibodies
Description	Rat IgG-Heavy and Light Chain Antibody
Target	Rat IgG-Heavy and Light Chain
Clonality	Polyclonal
Species/Host	Goat
Conjugation	Cy5
Reactivity	Rat
Predicted Reactivity	Bovine, Equine, Gallus, Human, Mouse, Rabbit, Sheep
Form/Appearance	Whole IgG
Concentration	0.5 mg/ml
Buffer/Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.09% Sodium Azide
Tested applications	FC, ICC, IF, IHC
Dilution range	IHC - 1:50 - 1:500; ICC - 1:50 - 1:500; IF - 1:50 - 1:500
Application notes	Format: Whole IgG
Antibody Type	Secondary Antibody
Storage	2 - 8°C

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Note

For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)