

## Product Datasheet

### Donkey anti-Rat IgG Heavy and Light Chain Antibody (HRP) (orb1531184)

|                             |                                                                                                     |
|-----------------------------|-----------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb1531184                                                                                          |
| <b>Category</b>             | Antibodies                                                                                          |
| <b>Description</b>          | Rat IgG-Heavy and Light Chain Antibody                                                              |
| <b>Target</b>               | Rat IgG-Heavy and Light Chain                                                                       |
| <b>Clonality</b>            | Polyclonal                                                                                          |
| <b>Species/Host</b>         | Donkey                                                                                              |
| <b>Conjugation</b>          | HRP                                                                                                 |
| <b>Reactivity</b>           | Rat                                                                                                 |
| <b>Form/Appearance</b>      | Whole IgG                                                                                           |
| <b>Concentration</b>        | 1 mg/ml                                                                                             |
| <b>Buffer/Preservatives</b> | Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.05% Pro-Clean 400                    |
| <b>Tested applications</b>  | ELISA, ICC, IHC, WB                                                                                 |
| <b>Dilution range</b>       | WB - 1:5,000 - 1:50,000; IHC - 1:200 - 1:5,000; ICC - 1:200 - 1:5,000; ELISA - 1:10,000 - 1:100,000 |
| <b>Application notes</b>    | Format: Whole IgG                                                                                   |
| <b>Antibody Type</b>        | Secondary Antibody                                                                                  |
| <b>Storage</b>              | 2 - 8°C                                                                                             |
| <b>Note</b>                 | For research use only                                                                               |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)