

Product Datasheet

ELISA Pro: Human IL-6 (orb1532209)

Catalog Number	orb1532209
Category	Tools
Description	The human IL-6 ELISA Pro kit provides all necessary reagents and solutions required to quantify human IL-6 in solution, e.g. in cell culture supernatant, and serum/plasma samples. The kit contains an ELISA plate pre-coated with monoclonal capture antibody, biotinylated detection antibody, Streptavidin-HRP, TMB substrate, ELISA standard and ELISA diluent. The ELISA diluent prevents false positive read-outs which may be caused by heterophilic antibodies found in serum and plasma samples. Presence of soluble IL-6 receptor does not significantly interfere with the assay.
Conjugation	HRP
Reactivity	Human, Primate
Form/Appearance	pre-coated plates
Tested applications	ELISA
Cross Reactivity	The monoclonal antibodies in this kit cross-react with IL-6 from non-human primates.
Range	10-3160 pg/ml
Storage	12 months
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558