

Product Datasheet

TACR2/NK2R Antibody (orb1533388)

Catalog Number	orb1533388
Category	Antibodies
Description	TACR2/NK2R Antibody
Target	TACR2 / NK2R
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human
Concentration	0.6 mg/ml
Buffer/Preservatives	PBS, 0.1% Sodium Azide
Purification	Immunoaffinity purified
Immunogen	Synthetic 20 amino acid peptide from C-Terminus of human Neurokinin A Receptor. Percent identity with other species by BLAST analysis: Human (100%); Gorilla, Gibbon (95%); Monkey, Marmoset (85%); Bat, Panda, Guinea pig (80%).
Tested applications	IHC, IHC-P
Dilution range	IHC, IHC-P (15 - 30 µg/ml)
Specificity	Human Neurokinin A Receptor. BLAST analysis of the peptide immunogen showed no homology with other human proteins.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Note For research use only

Expiration Date 12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)