

Product Datasheet

NPSR1/NPSR/GPR154 Antibody (orb1533913)

Catalog Number	orb1533913
Description	NPSR1/NPSR/GPR154 Antibody
Species/Host	Rabbit
Reactivity	Hamster, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 18 amino acid peptide from 2nd extracellular domain of human NPSR1. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Monkey, Marmoset, Mouse, Rat, Hamster (100%); Elephant, Panda, Bovine, Bat, Pig, Turkey, Chicken (94%); Dog, Opossum, Lizard (89%); Horse, Platypus (83%).
Target	NPSR1 / NPSR / GPR154
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Further information: Proteinase K treatment is required for antigen retrieval.
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (5 µg/ml)
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)