

Product Datasheet

GPR183/EBI2 Antibody (orb1534186)

Catalog Number	orb1534186
Category	Antibodies
Description	GPR183/EBI2 Antibody
Target	GPR183 / EBI2
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human, Mouse, Rat
Concentration	1 mg/ml
Buffer/Preservatives	PBS, 0.1% Sodium Azide
Purification	Immunoaffinity purified
Immunogen	Synthetic 20 amino acid peptide from 3rd cytoplasmic domain of human GPR183 / EBI2. Percent identity with other species by BLAST analysis: Human, Gorilla, Monkey, Mouse, Rat (100%); Marmoset, Elephant, Panda, Dog, Bovine, Horse, Rabbit, Pig, Opossum (95%); Platypus (90%); Hamster, Bat, Chicken (85%).
Tested applications	IHC, IHC-P
Dilution range	IHC, IHC-P (5 - 10 µg/ml)
Specificity	Human GPR183 / EBI2. BLAST analysis of the peptide immunogen showed no homology with other human proteins.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.



Brain, Medulla

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)