

Product Datasheet

FKSG89/MBOAT4 Antibody (orb1535022)

Catalog Number	orb1535022
Category	Antibodies
Description	FKSG89/MBOAT4 Antibody
Target	FKSG89 / MBOAT4
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human
Concentration	1 mg/ml
Buffer/Preservatives	PBS, 0.1% Sodium Azide
Purification	Immunoaffinity purified
Immunogen	Synthetic 14 amino acid peptide from internal region of human MBOAT4. Percent identity with other species by BLAST analysis: Human, Chimpanzee, Gorilla, Orangutan (100%); Gibbon, Monkey, Galago, Panda, Rabbit (93%); Mouse, Bat, Horse (86%).
Tested applications	IHC, IHC-P
Dilution range	IHC, IHC-P (5 - 10 µg/ml)
Specificity	Human MBOAT4. BLAST analysis of the peptide immunogen showed no homology with other human proteins, except METTL10 (57%).

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)