

Product Datasheet

GFPT1/GFAT Antibody (orb1535025)

Description	GFPT1/GFAT Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Equine, Gallus, Guinea pig, Human, Mammal, Monkey, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 12 amino acid peptide from N-Terminus of human GFPT1. Percent identity with other species by BLAST analysis: Human, Chimpanzee, Gorilla, Orangutan, Gibbon, Monkey, Galago, Marmoset, Mouse, Rat, Elephant, Panda, Bovine, Dog, Bat, Horse, Pig, Opossum, Guinea pig, Zebra finch, Chicken (100%); Turkey, Platypus (92%).
Target	GFPT1 / GFAT
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (2.5 - 5 µg/ml)
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)