

Product Datasheet

AXIN2/Axin 2 Antibody (orb1535572)

Catalog Number	orb1535572
Category	Antibodies
Description	AXIN2/Axin 2 Antibody
Target	AXIN2 / Axin 2
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human, Mouse
Concentration	1 mg/ml
Buffer/Preservatives	1X PBS, pH 7, 0.01% Thimerosal, 20% Glycerol. 1X PBS, pH 7, 0.01% Thimerosal, 20% Glycerol
Purification	Immunoaffinity purified
Immunogen	Recombinant protein encompassing a sequence within the center region of human Axin 2. The exact sequence is proprietary.
Tested applications	IHC, IHC-P, WB
Dilution range	IHC, IHC-P (1:100 - 1:1000), WB
Specificity	Human AXIN2 / Axin 2

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only



Anti-axin 2 antibody [N2C2], Internal used in IHC (Formalin-fixed paraffin-embedded sections).

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)