

Product Datasheet

STEAP1/STEAP Antibody (orb1535664)

Catalog Number	orb1535664
Category	Antibodies
Description	STEAP1/STEAP Antibody
Target	STEAP1 / STEAP
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human
Concentration	1.5 mg/ml
Buffer/Preservatives	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Purification	Caprylic acid and ammonium sulfate precipitation
Immunogen	Recombinant human Metalloreductase STEAP1 protein.
Tested applications	ELISA, IHC, IHC-P, WB
Dilution range	ELISA, IHC, IHC-P (5 µg/ml), WB
Specificity	Human STEAP1 / STEAP
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Note

For research use only



Immunofluorescent analysis of A431 cells using STEAP1 Antibody at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG (H + L)



Western blot All lanes: STEAP1 antibody at 1 µg/ml Lane 1: U87 whole cell lysate Lane 2: A431 whole cell lysate Lane 3: A549 whole cell lysate Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558