

Product Datasheet

HADHA Antibody (orb1536829)

Catalog Number	orb1536829
Category	Antibodies
Description	HADHA Antibody
Target	HADHA
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human, Mouse
Buffer/Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Purification	Affinity purified
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 484-763 of human HADHA (NP_000173.2). IAAVSKRPEKVIGMHYFSPVDKMQLLIITTEKTSKDTASAVAVGLKQGKVIIVVKDGPGFYT TRCLAPMMSEVIRILQEGVDPKKLDSLTSFGFPVGAATLVDEVGVDVAKHVAEDLGKVFG ERFGGGNPELLTQMVSKGFLGRKSGKGFYIYQEGVKRKDLNSDMDSILASLKLPPKSEVSS DEDIQFRLVTRFVNEAVMCLQEGILATPAEGDIGAVFGLGFPPCLGGPFRFVDLYGAQKIVD RLKKYEAAYGKQFTPCQLLADHANS PNKKFYQ
Tested applications	IF, IHC, IHC-P, WB
Dilution range	IF (1:50 - 1:200), IHC, IHC-P (1:100), WB (1:500 - 1:2000)

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom


Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC


68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Application notes	Further information: The predicted MW is 28kDa/83kDa, while the observed MW by Western blot was 82kDa.
Specificity	Human HADHA
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only



Immunofluorescence analysis of U2OS cells.



Western blot analysis of extracts of various cell lines.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558