

Product Datasheet

CRBN/Cereblon Antibody (orb1537041)

Catalog Number	orb1537041
Category	Antibodies
Description	CRBN/Cereblon Antibody
Target	CRBN / Cereblon
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human, Mouse, Rat
Concentration	1 mg/ml
Buffer/Preservatives	0.42% Potassium Phosphate, pH 7.3, 0.87% NaCl, 0.01% Sodium Azide, 30% Glycerol
Purification	Immunoaffinity chromatography
Immunogen	Recombinant full length protein of human CRBN
Tested applications	IHC-P, WB
Dilution range	IHC-P (1:100), WB (1:500 - 1:2000)
Application notes	Further information: Western blot analysis of CRBN expression in A549 (A), mouse spleen (B), rat brain (C) whole cell lysates.
Specificity	Recognizes endogenous levels of CRBN protein.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage


Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Expiration Date

12 months from date of receipt.



WB analysis of orb1537041.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)