

## Product Datasheet

### CBX3/HP1 Gamma Antibody (orb1537251)

<b>Catalog Number</b>	orb1537251
<b>Category</b>	Antibodies
<b>Description</b>	CBX3/HP1 Gamma Antibody
<b>Target</b>	CBX3 / HP1 Gamma
<b>Clonality</b>	Monoclonal
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Reactivity</b>	Human, Mouse
<b>Buffer/Preservatives</b>	PBS, pH 7.4, 0.02% Sodium Azide, 50% Glycerol
<b>Purification</b>	Purified
<b>Immunogen</b>	Recombinant protein
<b>Tested applications</b>	ICC, IHC-P, IP, WB
<b>Dilution range</b>	ICC (1:200), IHC-P (1:50), IP (1:500), WB (1:1000)
<b>Specificity</b>	Human CBX3 / HP1 Gamma
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)