

Product Datasheet

TGFB2/TGF Beta2 Antibody (orb1537317)

Catalog Number	orb1537317
Category	Antibodies
Description	TGFB2/TGF Beta2 Antibody
Target	TGFB2 / TGF Beta2
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human, Mouse
Buffer/Preservatives	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Purification	Caprylic acid and ammonium sulfate precipitation
Immunogen	Recombinant human Transforming growth factor beta-2 protein (aa304-413).
Tested applications	ELISA, IHC-P, WB
Dilution range	ELISA, IHC-P (10 µg/ml), WB
Specificity	Human TGFB2 / TGF Beta2
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Expiration Date 12 months from date of receipt.



Human Placenta: Formalin-Fixed, Paraffin-Embedded (FFPE)



Western blot All lanes: TGFB2 antibody at 3 μ g/ml Lane 1: Mouse heart tissue Lane 2: Mouse skeletal muscle tissue Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution Predicted band size: 48, 51 kDa Observed band size: 51 kDa

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558