

Product Datasheet

PTGS1/COX-1 Antibody (orb1537511)

Catalog Number	orb1537511
Category	Antibodies
Description	PTGS1/COX-1 Antibody
Target	PTGS1 / COX-1
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human, Mouse, Rat
Concentration	1 mg/ml
Buffer/Preservatives	0.42% Potassium Phosphate, pH 7.3, 0.87% NaCl, 0.01% Sodium Azide, 30% Glycerol. 0.42% Potassium Phosphate, pH 7.3, 0.87% NaCl, 0.01% Sodium Azide, 30% Glycerol
Purification	This antibody is purified by immunogen affinity chromatography
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-Terminus of human PGHS-1.
Tested applications	IHC-P, WB
Dilution range	IHC-P (1:100), WB (1:500 - 1:1000)
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



Human Skin: Formalin-Fixed, Paraffin-Embedded (FFPE)



Western blot analysis of PGHS-1 expression in HeLa (A); Raw264.7 (B); rat kidney (C) whole cell lysates.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)