

Product Datasheet

MANF/ARMET Antibody (orb1537855)

Catalog Number	orb1537855
Category	Antibodies
Description	MANF/ARMET Antibody
Target	MANF / ARMET
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human, Mouse, Rat
Concentration	1 mg/ml
Buffer/Preservatives	PBS containing 0.02% Sodium Azide
Purification	Peptide Column Affinity Chromatography
Immunogen	A 12 AA synthetic peptide near the C-Terminus of human MANF. The immunogen is located within the last 50 AA of MANF.
Tested applications	ELISA, IF, IHC, IHC-P, WB
Dilution range	ELISA, IF, IHC, IHC-P (5 µg/ml), WB (1 - 2 µg/ml)

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Application notes

Further information: MANF antibody can be used for detection of MANF by Western blot at 1 - 2 ug/mL. Antibody can also be used for immunohistochemistry starting at 2.5 ug/mL. For immunofluorescence start at 20 ug/mL. Antibody validated: Western Blot in mouse samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested. Western Blot: 20 kDa

Specificity

This antibody does not cross-react with CDFN.

Storage


Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note


For research use only

Expiration Date

12 months from date of receipt.



Immunofluorescence of MANF in Human Brain cells with MANF antibody at 20 ug/ml.



Western blot of MANF in mouse brain tissue lysate with MANF antibody at (A) 1 and (B) 2 ug/ml.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558