

Product Datasheet

DNMT3B Antibody (orb1538616)

Catalog Number	orb1538616
Category	Antibodies
Description	DNMT3B Antibody
Target	DNMT3B
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human
Concentration	1 mg/ml
Buffer/Preservatives	PBS (without Mg ²⁺ , Ca ²⁺), pH 7.4, 150 mM NaCl, 0.02% Sodium Azide, 50% Glycerol. PBS (without Mg ²⁺ , Ca ²⁺), pH 7.4, 150 mM NaCl, 0.02% Sodium Azide, 50% Glycerol
Purification	Immunoaffinity purified
Immunogen	Synthetic peptide from human DNMT3B. (AA Range: 1-50)
Tested applications	ELISA, IF, IHC, IHC-P, WB
Dilution range	IF (1:100 - 1:500), IHC, IHC-P (1:100), Peptide-ELISA (1:20000), WB (1:500 - 1:1000)
Application notes	Further information: Western Blot: Predicted MW 95 kDa

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Storage


Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note


For research use only



Anti-DNMT3B antibody IHC of human testis. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody dilution 1:100.



Immunofluorescence analysis of HeLa cells, using DNMT3B Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from HT-29 cells, using DNMT3B Antibody. The lane on the right is blocked with the synthesized peptide.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558