

Product Datasheet

LSM2/SnRNP Antibody (orb1539262)

Catalog Number	orb1539262
Category	Antibodies
Description	LSM2/SnRNP Antibody
Target	LSM2 / SnRNP
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG2a,k
Conjugation	Unconjugated
Reactivity	Human
Concentration	1 mg/ml
Buffer/Preservatives	PBS, pH 7.4, 0.1% Sodium Azide
Purification	Protein G purified
Immunogen	Recombinant human LSm2 (1-95aa) purified from E. coli.
Tested applications	ELISA, IHC, IHC-P, WB
Dilution range	ELISA, IHC, IHC-P (10 µg/ml), WB (1:1000 - 1:3000)
Specificity	Human LSM2 / SnRNP
Clone Number	AT2B2

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.



Anti-LSM2 / SnRNP antibody IHC staining of human testis. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.



Western Blot: The cell lysates (30 ug) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human LSm2 antibody (1:1, 000).

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558