

Product Datasheet

GAPDH-Loading Control mouse Monoclonal Antibody HRP Conjugated (orb1539606)

Catalog Number	orb1539606
Category	Antibodies
Description	GAPDH-Loading Control mouse Monoclonal Antibody HRP Conjugated
Clonality	Monoclonal
Species/Host	Mouse
Conjugation	HRP
Reactivity	Human, Mouse, Rat
UniProt ID	P04406
Tested applications	WB
Dilution range	WB=1:10000-1:50000
Antibody Type	Primary Antibody
Source	mouse
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)