

Product Datasheet

HSD17B3 Polyclonal Antibody PE Conjugated (orb1540076)

Catalog Number	orb1540076
Category	Antibodies
Description	HSD17B3 Polyclonal Antibody PE Conjugated
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	PE
Reactivity	Human, Mouse, Rat
UniProt ID	P37058
Tested applications	FC, IF
Dilution range	Flow-Cyt=1:50-200 IF=1:50-200
Antibody Type	Primary Antibody
Source	Rabbit
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)