

Product Datasheet

Beta Catenin (Phospho-Ser675) Antibody FITC Conjugated (orb1542438)

Catalog Number	orb1542438
Category	Antibodies
Description	Beta Catenin (Phospho-Ser675) Antibody FITC Conjugated
Clonality	Recombinant
Species/Host	Rabbit
Conjugation	FITC
Reactivity	Human, Mouse, Rat
UniProt ID	P35222
Tested applications	IF
Dilution range	Immunofluorescence: 1:50-200
Antibody Type	Primary Antibody
Source	Rabbit
Storage	Store at 4°C for 12 months.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)