

Product Datasheet

Integrin Alpha V + Beta 3 Antibody Cy5 Conjugated (orb1543035)

Description	Integrin Alpha V + Beta 3 Antibody Cy5 Conjugated
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Cy5
Tested Applications	IF
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Antibody Type	Primary Antibody
Source	Rabbit
Uniprot ID	P05106
Dilution Range	FCM:1:20-100IF:1:50-200
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)