

Product Datasheet

DNA PKcs ANTIBODY FITC Conjugated (orb1543307)

Catalog Number	orb1543307
Category	Antibodies
Description	DNA PKcs ANTIBODY FITC Conjugated
Clonality	Recombinant
Species/Host	Rabbit
Conjugation	FITC
Reactivity	Human, Mouse, Rat
UniProt ID	P78527
Tested applications	IF
Dilution range	Immunofluorescence: 1:50-200
Antibody Type	Primary Antibody
Source	Rabbit
Storage	Store at 4°C for 12 months.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)